Product Innovations

Dental - Labor | Dental Laboratory | Laboratoire Dentaire | Laboratorio Dentale | Laboratorio Dental
Zahnarzt - Praxis | Dental Practice | Cabinet Dentaire | Pratica Dentale | Práctica Dental

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All rights reserved.

Reproduction, also by extract and reproducing photos are only permitted with written authorization of EDENTA AG, Switzerland. All offers, orders and deliveries are subject to Edenta's „General sales and Delivery Terms“. We reserve the right to modify our range of products and their design as well as to deviate from the illustrations and data shown.

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Inhaltsverzeichnis | Index | Sommaire | Indice | Índice

	Page
DiaCut TURBO Diamonds (for quick preparation)	4
Occlusal shaper 972 (Diamond for all restoration)	5
K-Diamonds (for preparing of zirconia abutments)	6
Diamonds 801FGXL (for in-depth preparations)	7
Endo Access Bur C269GKS (for enlarging the pulp cavity)	8
Crown Cutter C35L (for separation of non precious and gold crowns)	9
Bone Cutter C254A (for preparation of bone and tooth structure)	10
Bone Cutter C141AC (enhanced with ZrN hard coating)	11
Perforated Diamond Strips PDS2 / PDS3 / PXDS2 (for contouring)	12
Initial Burs 186RF	13
Pilot Burs 210L	13
Trepan Burs 229, 229L	14
Cutters with internal cooling 81IK, 167IK, 210IK,	15
StarTec Polisher RA (2-step Polishing System for Lithium-Disilicat)	16
Prophylaxis paste	17
StarTec Polisher HP (2-step Polishing System for Lithium-Disilicat)	18
Preparing of unsintered zirconia restorations	19
TC Cutter VZ65 (for working PMMA)	20
ExaStar Polisher HP (2-step Polishing System for PMMA)	21
TC Cutter VZ86 (for processing plaster)	22
VOLCANO TC Cutter (enhanced with DLC hard coating)	23

DiaCut TURBO Diamonds

VT830L

Birne, lang
Long Pear
Poire longue

⌚ max. 450.000 / opt. 70.000-220.000

📦 5

Shank	L mm	ISO	Order No.	4,0	5,0
Size				012	014
∅ ₁₀ mm					
FG	• VT806 314 239 534...		VT830L.314...	012	014

VT850

Konus, rund
Round End Taper
Cone, bout arrondi

⌚ max. 300.000 / opt. 70.000-140.000

📦 5

Shank	L mm	ISO	Order No.	10,0
Size				016
∅ ₁₀ mm				
FG	• VT806 314 199 534...		VT850.314...	016

VT879

Torpedo
Beveled Cylinder
Torpille

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DiaCut TURBO Diamonds

High-quality diamond instruments with a high cutting capacity for atraumatic, quick preparation. The combination of diamond-coated spirals and multilayer diamond matrix produces a smooth preparation surface without grooves. The non-diamond coated spiral section guarantees reliable cooling to avoid pulp damage due to excessive heat build-up. Large chip escape channels ensure rapid removal of the debris and consequently optimum self-cleaning of the instruments. The special TiN coating retains the diamonds tightly, combined with a high instrument service life. A specially adapted production technology and high-grade raw materials guarantee constant product quality.

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
Occlusal shaper 972

A central aspect of restorative dentistry is the design of the occlusal surfaces. These are not only important from a masticatory function point of view but also for the stability of the entire stomatognathic system. The occlusal surfaces of a denture tooth always orientate according to the natural shape of the tooth with functional cusps and fissures, which in turn must integrate with those of the adjacent and opposing teeth with regard to the anatomical shape. The Occlusal shaper 972, which was specially developed for occlusal contouring, is suitable for all filling materials and also all-ceramic ZrO₂ material for reliable contouring of the cusps and fissures using only one instrument. The round shape facilitates contouring of cusps and fissures. The curvature prevents cutting of notches. The fine and superfine grit produce a perfect surface for the subsequent polish.

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
K879

Torpedo
Modified Beveled Cylinder
Torpille

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K-Diamonds with 80µm special grit

Zircon oxide abutments are the perfect solution for all-ceramic superstructures. The abutments are contoured using the newly developed K-Diamond instruments. The 80 µm special grit prevents damage of the zircon oxide all-ceramic. Crack formation and breaking out of crystals from the zircon oxide is prevented. The high stability of the special grit produces an outstanding service life of the K-Diamond instruments. The new electroplating bonding prevents the diamond grit breaking out prematurely. This creates an unrivalled abrasive capacity, perfectly, quickly and reliably.

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
Diamonds 801FGXL

Diamond 801 with an extra-long shank FGXL and different grit sizes. During endodontic procedures and with cavity preparations the diamond instruments with an extra-long shank facilitate the necessary preparations in deeper areas. The FGXL shank improves the visibility in the deep areas of the cavity. The different grit sizes enable individual preparation.

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
Endo Access Bur C269GKS

Special endodontic bur with a long, conical blade configuration and cross-cut used for preparing a funnel-shaped access in the pulp chamber. The cross-cut guarantees quiet, efficient removal of tooth structure. The rounded, non-cutting bur tip (safe end) prevents penetration into the root canal walls and pulp chamber floor. Application after removal of the pulp chamber roof using the endodontic bur 802L for extending the access opening to ensure successful endodontic treatment.

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
Endo Access Bur 802L for removal of the pulp chamber roof and extending the access opening.

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Crown Cutter C35L for metal crowns

Specially coordinated for the separation of non-precious metal and gold crowns with an improved cutting capacity and service life compared with conventional crown cutters. The special compact blade configuration provides the C35L with controlled cutting behaviour and high cutting efficiency due to the optimally adapted tip-transversing blade with cross-cut. Use in the red contra-angle at 160,000 rpm; can also be used in a turbine.

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

Drehzahl opt. 120'000 – 160'000 upm, Anpresskraft: 0.5 N

Recommended use:

- Cut open **ceramic veneer** using diamond instrument.
- Separate **metal framework** by applying several times at less than 45° and machining small pieces.

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


Bone Cutter C254A

Bone cutter used for preparing bone and tooth structure. The special A-cut blades guarantee very quiet instrument running at the highest cutting capacity with maximum conservation of the tooth structure to be prepared. The long, slimline neck provides excellent visibility during preparation.

Indications: Exposing and removing impacted wisdom teeth, apicectomies, preparation of bone structure and a bone lid, separating teeth and tooth roots.

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
Bone Cutter C141AC

Bone cutter used for gentle preparation of the bone structure. The new optimised A blades guarantee very quiet instrument running with the highest cutting efficiency and minimum trauma of the bone structure to be machined.

The cutter meets the high demands placed on a high-quality surface during surgical preparation of the bone, due to the ZrN (zircon nitride) hard coating. The transitional blade also enables axial insertion in the bone structure. The newly developed blade design guarantees controlled preparation and a smoother bone structure surface for good wound healing.

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Perforated Diamond Finishing Strips

Diamant Finier Streifen mit Sägeverzahnung und Perforation, einseitig belegt mit Diamantkörnung, sterilisierbar.
 Diamond finishing strips with serrated edge and perforation, coated on one side with diamond grit, sterilisable.
 Strips à finir diamantés perforés à profil denté (scie) recouverts de grains de diamant sur une face, stérilisables.
 Strisce diamantate per finitura con dentatura a sega, con diamantatura unilaterale, sterilizzabili.
 Tiras de acabado diamantadas con dentado de sierra y perforadas, con cobertura unilateral de granulado de diamante, esterilizables.

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Perforated Diamond Strips are designed for complete control during interproximal reduction, shaping, and contouring. Strips allow easy access and precise manual enamel reduction resulting in a smooth, natural finish. Perforated design assists in debris removal, provides improved visibility, control, and flexibility. Single-sided stain-less steel strip resistant to breaking and stretching, reusable, autoclavable.

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

Initial Bur 186RF

The initial bur is used as the first instrument during surgical drilling procedures to place a drill mark and thus prevent the subsequent spiral drill skidding, in particular with tapering alveolar ridges.

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

Pilot Bur 210L

The efficiently cutting pilot bur is suitable as a pre-drill for all commonly used implant systems for initial preparation of the implant site. The specially designed bur tip and the good axial drilling properties enable precise drilling. The large chip escape channels guarantee rapid chip removal. The lasered depth marks at intervals of 2 millimetres provide exact visualisation of the penetration depth.

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

Trepan bur 229, 229L

Specially suitable for use in oral surgery and in implantology for accurate, easy preparation of cylindrical bone grafts and bone extraction (bone lid during apicectomies). The good axial drilling properties enable precise drilling. The efficiently cutting bone blades allow atraumatic treatment. The lasered depth markings 2 millimetres apart provide exact visualisation of the penetration depth.

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

Bone cutter with internal cooling (IK) for gentle preparation of bone structure, stainless steel.

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Pilot bur with internal cooling (IK) for implantology, stainless steel with depth marking = 2 mm

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Bone cutter and drills with internal cooling

In surgical work bone preparation must be completed with as little heat build-up as possible to ensure reliable healing. The special drilling and cutting instruments with internal cooling are suitable for efficient cooling of the bone.

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StarTec Polisher RA

Specially developed for gentle polishing of lithium disilicate glass ceramic. After milling in a CAD/CAM machine ceramic restorations have a rough surface. To reduce abrasion of the opposing dentition the surface roughness must be smoothed. The StarTec polishers guarantee finishing and polishing of the ceramic surface and therefore minimal surface roughness.

Advantages of StarTec polishers: 2-step diamond polishing system for reliable, gentle polishing of glass ceramics, silicate ceramics, feldspar ceramics, lithium disilicate ceramics and leucite ceramics. Optimal finishing and polishing in only 2 steps. A newly developed, highly condensed diamond grit binder produces an improved surface quality. The specially developed polisher matrix guarantees an optimum surface sheen without the use of polishing paste.

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

Prophylaxis paste

Green, medium. RDA 170, corresponds to 70 microns.
Contents 200 x 2g pot and 1 ring holder. Use for severe discoloration and plaque accumulation.

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Prophylaxis paste for maximum cleaning of the tooth surface

Healthy teeth and healthy gingiva – with the correct prophylaxis. The new prophylaxis paste with 1.23% fluoride and optimum cleaning power and abrasion for removing tooth accretion and discoloration of the dental enamel to prevent renewed fixation of bacteria.

Also used for polishing following calculus removal, removal of orthodontic brackets as well as the cleaning of filling materials. Specially formulated paste to avoid splashing during polishing. Pleasant peppermint taste. Easy dispensing and handling thanks to hygienically packed, easy-to-use single doses 2g pot with ring holder.

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StarTec Polisher HP

2-step diamond polishing system for reliable and gentle preparation of zirconium dioxide, lithium disilicate, zircon-reinforced lithium silicate (ZLS), hybrid-ceramic and all commonly used veneering porcelains. After milling in a CAD/CAM machine, ceramic restorations have a surface roughness. To reduce the abrasion on the opposing tooth the surface roughness must be reduced. The StarTec polishers guarantee finishing and polishing of the ceramic surface and therefore minimum surface roughness.

Advantages of StarTec polishers: 2-step diamond polishing system for reliable and gentle polishing of glass-ceramics, silicate ceramics, feldspar ceramics, lithium disilicate and leucite ceramics. Optimum finishing and polishing in only 2 steps. The newly developed, highly compressed binding of the diamond grit produces an improved surface quality. The specially developed polisher matrix guarantees an optimum surface polish without polishing paste and makes a renewed glaze firing unnecessary after adjustment.

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Trimming unsintered zirconia copings before sintering/stain sintering

A prerequisite for a high service life of all-ceramic restorations is material-specific and therefore material-friendly preparation of the unsintered ZrO₂ restoration. Optimum preparation of “soft” zircon oxide (ZrO₂ partially sintered copings) is completed as described in the following (1/2/3).

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TC Cutter, 0175.023HP Separation of the retention bars. Optimum speed 20.000 rpm

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TC Cutter, 0830.023HP – Trimming attachment areas of the retention bars and adjustment of uneven areas. Optimum speed 20.000 rpm

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Polisher, 0664HP-0669HP – Smoothing of the surfaces to ensure colouring liquids penetrate to an optimum depth. Optimum speed 10 – 15.000 rpm

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Tungsten carbide cutter (VZ65) for trimming PMMA materials

Specially developed tungsten carbide cutters for finishing and reworking PMMA restorations. With a new blade geometry (VZ65) for gentle trimming of thermoplastic PMMA materials, which have the characteristic of smearing quickly. The new type of blade enables quick material removal with a sharp cut without heating the material, does not snag and produces smooth material surfaces, allowing them to be easily polished.

Applications: For separating retention bars (0165.023HP), trimming the attachment areas of the retention bars and for adjustment of uneven areas in the shape. **Optimum speed** 20.000 rpm, only apply light pressure during preparation.

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
ExaStar Polisher for polishing PMMA

ExaStar						
L mm	3,0	3,0	15,5	3,0	3,0	15,5
Size \varnothing 1/10 mm	170	220	055	170	220	055
\varnothing opt.	12.000	12.000	12.000	6.000	6.000	6.000
\varnothing max.	15.000	15.000	15.000	15.000	15.000	15.000
Order No.	0612HP	0610HP	0614HP	0622HP	0620HP	0624HP
ISO No. 652 104...	372 524 170	303 524 220	243 524 055	372 514 170	303 514 220	243 514 055
Stufe / Step / Etape	1			2		

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ExaStar polisher for polishing PMMA materials

Specially developed 2-step polishing system with hybrid grit and matrix binder for the reliable and gentle polishing of thermoplastic PMMA materials. Careful polishing reduces plaque accumulation and is thus a prerequisite for good aesthetics. After milling the PMMA restoration in the CAD/CAM machine, the surface roughness must be reduced. ExaStar polishers guarantee a finish and glaze polish for an optimum surface lustre. The new polisher matrix in combination with the hybrid grit enables an optimum surface polish without the use of polishing paste. Polishing easily and quickly produces an aesthetic result, so that the PMMA restoration can be fitted immediately after polishing.

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



Tungsten Carbide Cutter (VZ86) for trimming plaster models

The newly developed VZ86 ensures reliable trimming of plaster. The phase cut guarantees high material removal and a smooth material surface. The anticlockwise blade configuration assures reliable trimming, as the cutter shank does not pull out of the chuck, even with high material removal. The large chip escape channels guarantee optimum chip transport and prevent clogging of the cutter. Thanks to the phase cut the cutter is particularly quiet running, is easy to manage and does not catch.

Optimum speed 10.000 - 15.000 rpm, only apply light pressure during trimming.

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
Volcano - Tungsten Carbide Cutter with DLC hardening for increasing the performance, combined with a high cutting capacity and maximum service life. The DLC wear protection creates a hard cutter with a smooth, porosity-free cutting surface. This increases the cutter hardness and reduces the frictional resistance. The wear resistance is increased and premature breaking out of the cutter blades is prevented, which considerably extends the service life.

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For the preparation of titanium, non-precious metal alloys, precious metals, model dental stone, CrCo, veneering resins, denture acrylics.

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For trimming titanium and titanium alloys.

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For fine finishing all alloys and composites.

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